

PINK beamline

PShell custom functions

Table of Contents

scan.....	2
continuous.....	2
filter_scan.....	2
gap_scan.....	3
line.....	3
spot.....	3
zigzag.....	4
sample_scan.....	4
pink.....	5
bl_snapshot_print.....	5
bl_snapshot_save.....	5
filter<n>.....	5
gap.....	5
ml_<energy>ev.....	5
rename_sample.....	6
shutter_hard_CLOSE.....	6
shutter_hard_OPEN.....	6
valveCLOSE.....	6
valveOPEN.....	6
AU1.....	6
eiger_set_background.....	6
plot.....	7
scatter.....	7
spec_atm_go2_stdby_position.....	7
bpm.....	8
BPM<n>_Horizontal_Profile.....	8
BPM<n>_Vertical_Profile.....	8
blade.....	9
diagnostic_chamber_blade_scan.....	9
Sample_Env_Blade_Scan.....	9
Slit_Scan.....	9

scan.

continuous

scan.continuous(detector, det_exposure, sample_exposure, X0, X1, dX, Y0, Y1, passes, sample, linedelay)

detector:	Detector selection
det_exposure:	Detector's exposure time
sample_exposure:	Sample's exposure time per pass
X0:	start location on X axis
X1:	End location on X axis
dX:	Distance between vertical lines
Y0:	Start location on Y axis
Y1:	End location on Y axis
passes:	Number of passes
sample:	Name of the sample
Linedelay:	Time delay between vertical lines




Fig. 1 - Continuous scan movement diagram (beam perspective looking downstream).

Performs a scan on a moving sample. The goal is to maximize detector's exposure time and minimize sample degradation. The speed of the sample motion is calculated using the defined "sample_exposure" time and a vertical beam size of 50 micrometers.

filter_scan

scan.filter_scan(filters, start, end, step, exposure)

Filters:	Filter selection
start:	Start position on filter
end:	End position on filter
Step:	Distance between point (float) or number of steps in interval (int)
Exposure:	Electrometer's reading time

Performs a filter position scan against Izero. The beam intensity readout comes from the sum of 4 quadrant diodes used for beam positioning monitoring.

gap_scan

scan.gap_scan(start, end, step, exposure, fit)

start:	starting undulator gap
end:	end undulator gap
step:	Step (float) or number of points to acquire (int)
exposure:	Electrometer acquisition time
fit:	If True, fits data to a gaussian curve with a exponential background

Performs an Undulator gap scan. The beam intensity readout comes from the sum of 4 quadrant diodes used for beam positioning monitoring.

line

scan.line(detector, exposure, Y0, dY, Ypoints, passes, sample)

Detector:	Detector selection
exposure:	detector's exposure time
Y0:	Start location on Y axis
dY:	distance between points
Ypoints:	number of points
passes:	number of passes
sample:	name of the sample

Performs multiple exposures on a sample along a vertical line.

spot

scan.spot(detector,exposure, images, sample)

detector:	detector selection
exposure:	detector's exposure time
images:	number of exposures
sample:	name of the sample

Performs multiple exposures on a single location on sample.

zigzag

scan.zigzag(detector, exposure, X0, dX, Xpoints, Y0, dY, Ypoints, passes, sample, linedelay)

detector:	detector selection
exposure:	detector's exposure time
X0:	Start location on X axis
dX:	distance between columns in the X axis
Xpoints:	number of columns in the X axis
Y0:	Start location on Y axis
dY:	distance between rows in the Y axis
Ypoints:	number of rows in the Y axis
passes:	number of passes
sample:	name of the sample
linedelay:	time delay between each line exposure

Performs multiple exposures with a grid pattern on a sample along X and Y lines.

sample_scan

scan.sample_scan(axis, detector, start, end, step, exposure)

axis:	axis selection: vertical or horizontal
detector:	detector selection
start:	Start position on sample
end:	End position on sample
step:	Distance between point (float) or number of steps in interval (int)
exposure:	detector's exposure time

Performs a sample position scan against the detector's peak-count-per-second value. Useful to find sample edge positions.

pink.

bl_snapshot_print

pink.bl_snapshot_print()

Prints on screen the position, temperature and pressure of relevant components for the Pink beamline.

bl_snapshot_save

pink.bl_snapshot_save()

Prints on screen the position, temperature and pressure of relevant components for the Pink beamline.

filter<n>

pink.filter1(pos), filter2(pos), filter3(pos)

pos:	position of the filter
-------------	------------------------

Moves filter to the requested position. All filters are located inside the diagnostic chamber.

gap

pink.gap(value)

value:	no value or undulator gap (mm)
---------------	--------------------------------

Prints the gap position if no value is given; Or sets the gap on the undulator.

ml_<energy>ev

pink.ml_2300ev(), ml_3000ev(), ml_4000ev(), ml_5000ev(), ml_6300ev(), ml_6800ev(), ml_7300ev(), ml_8000ev(), ml_9500ev()

Sets up the correct stripe on the multilayer mirror based on the desired beam energy.

rename_sample

pink.rename_sample(path, newstring)

path:	variable path (can be retrieved by right clicking on the variable using the data explorer tab)
newstring:	new string

Rename the sample on a acquire data file (.h5 data file)

shutter_hard_CLOSE

shutter_hard_OPEN

pink.shutter_hard_OPEN(), shutter_hard_CLOSE()
Open/close hard beam shutter.

valveCLOSE

valveOPEN

pink.valveOPEN(vnum), valveCLOSE(vnum)

vnum:	Valve number according to diagram on control panel
--------------	--

Open/close the vacuum valve.

AU1

pink.AU1(top, bottom, wall, ring):

top,bottom,wall,ring:	Aperture defined in mm.
------------------------------	-------------------------

Set the specified aperture on AU1.

eiger_set_background

pink.eiger_set_background():

:	no argument
----------	-------------

Execute sequence to set the correct background information on eiger detector after reboot of detector's IOC.

plot

pink.plot(filepath)

filepath:	Copied file location from the data tab
------------------	--

Plots relevant data from h5 file from the data tab. To copy the filepath, click first on the desired file on the left window, then right-click the name of the file on the middle window. Click on "Copy link to clipboard", then past it as the filepath.

scatter

pink.scatter(pos)

pos:	position of the scatter filter
-------------	--------------------------------

Moves scatter filter to the requested position. The scatter filter is located inside the diagnostic chamber.

spec_atm_go2_stdby_position

pink.spec_atm_go2_stdby_position()

:	no argument
---	-------------

Moves ATM spectrometer to the standby position for detector replacement.

bpm.

BPM<n>_Horizontal_Profile

BPM<n>_Vertical_Profile

bpm.BPM<n>_Vertical_Profile()

Provides information about the beam dimensions on the selected BPM (Beam Position Monitor)

blade.

diagnostic_chamber_blade_scan

blade.Diagnostic_Chamber_Blade_Scan(start, end, step, exposure)

start:	Start position on filter's blade
end:	End position on filter's blade
step:	Distance between points (float) or number of steps in interval (int)
exposure:	detector's exposure time

Perform a blade scan with blade inside the diagnostic chamber. After scan, the script attempts to perform a gaussian fit to derivative of the acquired data.

Sample_Env_Blade_Scan

blade.Sample_Env_Blade_Scan(start, end, step, exposure)

start:	Start position on sample holder's blade
end:	End position on sample holder's blade
step:	Distance between points (float) or number of steps in interval (int)
exposure:	detector's exposure time

Perform a blade scan with blade on the sample holder. After scan, the script attempts to perform a gaussian fit to derivative of the acquired data.

Slit_Scan

blade.Slit_Scan(slit, start, end, step, exposure)

slit:	slit selection
start:	Start position on sample holder's blade
end:	End position on sample holder's blade
step:	Distance between points (float) or number of steps in interval (int)
exposure:	detector's exposure time

Perform a blade scan with the selected slit. After scan, the script attempts to perform a gaussian fit to derivative of the acquired data.